


# IChem: A Toolkit for detecting, comparing and predicting protein-ligand interactions



**Jérémy DESAPHY, Franck DA SILVA, Guillaume BRET and Didier ROGNAN**

*Laboratoire d'Innovation Thérapeutique, UMR 7200 CNRS-Université de Strasbourg, F-67400 ILLKIRCH*

Email: rognan@unistra.fr



## Literature Corner

Please have a look at these articles for detailed information on the basic principles and concepts underlying IChem usage.

Marcou, G. and Rognan, D. (2007) Optimizing fragment and scaffold docking by use of molecular interaction fingerprints. <i>J. Chem. Inf. Model.</i> , <b>47</b> , 195-207	
Desaphy, J., Azdimousa, K., and Rognan, D. (2012) Comparison and druggability prediction of protein-ligand binding pockets from pharmacophore-annotated shapes., <i>J. Chem. Inf. Model.</i> , <b>52</b> , 2287-2299	
Desaphy, J., Ducrot, P., Raimbaud, E. and Rognan, D. (2013) Encoding protein-ligand interaction patterns in fingerprints and graphs. <i>J. Chem. Inf. Model.</i> , <b>53</b> , 623-637	
Desaphy, J. and Rognan, D. (2014) scPDBFrag: a database of protein-ligand interaction patterns for bioisosteric replacements. <i>J. Chem. Inf. Model.</i> , <b>54</b> , 1908-1918	
Gabel, J., Desaphy, J. and Rognan, D. (2014) Beware of machine learning-based scoring functions - On the danger of developing black boxes. <i>J. Chem. Inf. Model.</i> , <b>54</b> , 2807-2815	
Da Silva, F., Desaphy, J., Bret, G. and Rognan, D. (2015) IChemPIC: A Random Forest Classifier of Biological and Crystallographic Protein-Protein Interfaces. <i>J. Chem. Inf. Model.</i> , <b>55</b> , 2005-2014	
Slynko, I., Da Silva, F., Bret, G. and Rognan, D. (2016) Docking pose selection by interaction pattern graph similarity: application to the D3R grand challenge 2015. <i>J. Comput.-Aided Mol. Des.</i> , <b>30</b> , 669-683.	
da Silva Figueiredo Celestino Gomes, P., Da Silva, F., Bret, G. and Rognan, D. (2018) Ranking docking poses by graph matching of protein-ligand interactions: lessons learned from the D3R Grand Challenge 2. <i>J. Comput.-Aided Mol. Des.</i> , <b>32</b> , 75-87.	
Da Silva, F., Desaphy, J. and Rognan D. (2018) IChem: A Versatile Toolkit for Detecting, Comparing, and Predicting Protein-Ligand Interactions. <i>ChemMedChem</i> , <b>13</b> , 507-510.	
Tran-Nguyen, V.K., Da Silva, F., Bret, G. and Rognan, D. (2019) All in One: Cavity Detection, Druggability Estimate, Cavity-Based Pharmacophore Perception and Virtual Screening. <i>J. Chem. Inf. Model.</i> , <b>59</b> , 1, 573-585	
Koensgen, F., Da Silva, F., Rognan, D. and Kellenberger, E. (2019) Unsupervised Classification of G-Protein Coupled Receptors and Their Conformational States Using IChem Intramolecular Interaction Patterns, <i>J. Chem. Inf. Model.</i> , <b>59</b> , 3611-3618	
Da Silva, F., Bret, G., Teixeira, L., Gonzalez, C.F. and Rognan, D. (2019) Exhaustive repertoire of druggable cavities at protein-protein interfaces of known three-dimensional structure. <i>J. Med. Chem.</i> , <b>2019</b> , 62, 9732-9742	
Tran-Nguyen, V.-K., Jacquemard, C. and Rognan D. (2020) LIT-PCBA: An unbiased dataset for machine learning and virtual screening. <i>J. Chem. Inf. Model.</i> , <b>60</b> , 4263-4273	
Tran-Nguyen, V.-K., Bret, G. and Rognan, D. (2021) Accuracy of fast scoring functions to predict high-throughput screening data from docking poses: The simpler the better. <i>J. Chem. Inf. Model.</i> , <b>67</b> , 2788-2797	
Volkov, M., Turk, A.J., Drizard, N., Martin, M., Hoffmann, B. Gaston-Mathé, Y. and Rognan, D. (2022) On the frustration to predict binding affinities from protein-ligand structures with deep neural networks. <i>J. Med. Chem.</i> , <b>65</b> , 7946-7958	



## What's new ?

IChem v 5.2.9 offers the possibility to automatically generate cavity-based pharmacophores from the sole 3D structure of a target of interest. This new functionality is embedded in the Volsite module and enables the user to define structure-based pharmacophores in three possible formats:

- .chm for use in Biovia DiscoveryStudio or PipelinePilot
- .pml for use in LigandScout
- .mol2 for use in Shaper2 (in-house ligand alignment tool)

For more details about the procedure, please have a look at:


Tran-Nguyen VK, Da Silva F, Bret G, Rognan D.  
J Chem Inf Model. 2019, 59, 1, 573-585

[All in One: Cavity Detection, Druggability Estimate, Cavity-Based Pharmacophore Perception and Virtual Screening.](#)



## Installation

IChem is provided as a zipped archive file (**IChem.tgz**) containing the following material:



- **lib**: a directory containing necessary template files
- **test**: a directory containing some test input/output files
- **Ichem.lic** a license file
- **IChem**: A 64-bit Linux executable (CentOS 7.1)
- **User\_Guide.pdf**: this manual

Source some environment variables and untar the IChem distribution in the IChem root directory

For csh users: *setenv ICHEM\_DIR /your\_IChem\_root\_directory*  
*setenv ICHEM\_LIC \$ICHEM\_DIR/IChem.lic*  
*setenv ICHEM\_LIB \$ICHEM\_DIR/lib*  
*cp IChem.tgz \$ICHEM\_DIR*  
*cd \$ICHEM\_DIR; tar -cvfz IChem.tgz*

for bash users: *ICHEM\_DIR="/your\_IChem\_root\_directory"; export ICHEM\_DIR*  
*ICHEM\_LIC= \${ICHEM\_DIR}/IChem.lic; export ICHEM\_LIC*  
*ICHEM\_LIB= \${ICHEM\_DIR}/lib; export ICHEM\_LIB*  
*cp IChem.tgz \${ICHEM\_DIR}*  
*cd \${ICHEM\_DIR}; tar -cvfz IChem.tgz*

### Note to users:

IChem commands with options/tool/arguments will be displayed in *italic* characters after the ">" prompt

Input/output filenames will be displayed in **bold violet** characters

Terminal output will be displayed with a **gray background**

A copy of all test input/output files is given in the test directory of the IChem distribution. Before using IChem, please take the time to read the description of the technology. Articles to read will be mentioned by the following icon:





## A short introduction to IChem

IChem is a multi-task program for detecting, analyzing and comparing protein-ligand interactions. It is composed of several tools:

Tool	Purpose
<b>Binding mode tools</b>	
<b>IFP</b>	Detect and save interactions as fingerprints (binding-site dependent)
<b>ints</b>	Detect and save interaction patterns( binding-site independent)
<b>grim</b>	Detect and save interactions as graphs
<b>Generic tools</b>	
<b>genkey</b>	License file generator (privileged usage)
<b>license</b>	Output license details
<b>pdbconv</b>	Convert PDB into MOL2 files
<b>realign</b>	Rotation/translation of atomic coordinates
<b>sims</b>	Fingerprint comparison
<b>utils</b>	Miscellaneous (buried surface area, ligand fragmentation)
<b>Cavity detection</b>	
<b>Volsite</b>	Cavity detection and druggability/ligandability prediction
<b>Protein-Protein interfaces</b>	
<b>DetectPPI</b>	Cavity detection and druggability/ligandability prediction



## IChem Usage

The correct syntax for using IChem is:

---

```
> IChem options tool arguments
```

---

Alternatively, options, tool and arguments can be stored in a parameter file (any name) that can be called from IChem with the `-F` option:

---

```
> IChem -F parameter_file
```

---

Example of a parameter\_file content

```
options_1 tool_1 arguments_1
options_2 tool_2 arguments_2
...
options_n tool_n arguments_n
```

Multiple IChem commands can be run from a single parameter file, just by adding for every line a novel list of options/tool/arguments

### Verbose and debugging mode:

IChem has been written in order to simplify standard output. If you want more details than that saved in the regular output, please use the `--verb` option:

---

```
> IChem --verb options tool arguments
```

---

Moreover, a debugging mode is available in IChem by calling the `--debug` option with one or more of the following values, separated by a comma:

```
READ  MOL2 file parser
DINT  Binding mode detection
GRID  Three-dimensional grid
MOLD  Molecular analysis
```

---

```
> IChem -debug READ,DINT,MOLD options tool arguments
```

---



## IChem menu

Just typing **\$ICHEM\_DIR/IChem** (or just **IChem** if \$ICHEM\_DIR is sourced in your path) gives you access to the full IChem menu

---

**> IChem**

---

### KEY Generator

**genKey** *year month day allowed\_tools*

year :                    Expiring year (4 digits : 2014)  
month:                    Expiring month (2 digits : 09)  
day :                      Expiring day (2 digits : 23)  
allowed\_tools: each value is separated by space and must have a value of either 1 or 0. Must follow this order  
1st:                        Licence Key generator  
2nd:                        Molecular realignment  
3rd:                        BSA Calculation  
4th:                        IFP generator  
5th:                        Interaction detection  
6th:                        Graph Interaction Matching  
7th:                        VolSite  
8th:                        PDB to MOL2 conversion  
9th:                        Patching MOL2 conversion  
10th:                      Utils  
11th:                      Sims  
12th:                      Scoring  
13th:                      Detection and analysis of PPI

Example :

```
genKey 2014 1 28 0 0 1 0 0 0 0 1 0 1 1 1
      |<-DATE->|<-----TOOLS----->
```

## realign - Molecular alignment

**realign** *rigidM mobilM applied1 applied2*

||-> rigidM : reference molecule to apply alignment to  
 ||-> mobilM : comparison molecule to apply alignment from  
 ||-> applied: molecule to apply rotation/translation to

### [General options]

-gmatch N (NAME) Use graph matching to align  
 NAME Atom Name matching  
 ATMN Atomic Name matching  
 MOL2 MOL2 Type matching  
 CALP CAlpha Atom matching (protein only)  
 --wMob Also outputs the aligned mobilM  
 -rule R

By default, the program will perform an atom by atom match, without taking care of what kind of atom it match. If you want to perform a match by regarding only some atoms, this index\_string is here to do so

ex : -i '2-3|1-6|23-160' Will match the second atom from the reference with the third from the comparison, the first with the sixth ...

## IFP - Interaction FingerPrint

**IFP** *protein ligand*

**IFP** *protein ligand ligand\_ref*

### [General options]

-name N (LIG) Name of the fingerprint -Default: Name of the ligand  
 --polar Detect and output only polar interactions  
 --metal Detect and output only metal interactions  
 --extended Include within the fingerprint:  
 |--> Metal/Acceptor interaction  
 |--> Weak Hydrogen bonds  
 |--> PI-Cation interactions

### [testing options]

-D\_Hb N (3.5) Hbond length (Angstroem)  
 -D\_Hyd N (4.5) Hydrophobic length (Angstroem)  
 -D\_Io N (4.0) Ionic length (Angstroem)  
 -D\_Me N (2.8) Metal/Acceptor length (Angstroem)  
 -D\_Ar N (4.0) Aromatic interaction length (Angstroem)  
 -D\_Pic N (4.0) Pi cation interaction length (Angstroem)  
 -a\_H N (Pi) HBond angle (rad)  
 -at\_H N (Pi/3) HBond tolerance angle (rad)  
 -a\_ArFF N (Pi) Aromatic Face to Face interaction angle (rad)  
 -at\_ArFF N (Pi/6) Aromatic Face to Face tolerance angle (rad)  
 -a\_ArEF N (Pi/2) Aromatic Edge to Face interaction angle (rad)  
 -at\_ArEF N (Pi/3) Aromatic Edge to Face tolerance angle (rad)  
 -a\_Pic N (Pi) Pi cation interaction angle (rad)  
 -at\_Pic N (Pi/6) Pi cation tolerance angle (rad)  
 --ligD Print all possible ligand interactions

Please note that ligand file can be a multi-mol2 file



## INTERACTION GENERATOR

### *ints prot lig out*

-type (CENT)	Alter positionning output, multiple values are allowed, separated by space
PROT	InterPROT positionning
LIG	InterLIG positionning
CENT	Centered positionning
MERG	Merged all 3 above
-fgps	Fingerprint format
STD	Standard (1 0 21 0 0 3)
SVM	SVM format (1:1 3:21 6:3)
CMP	Compressed (1 [1 21 [2 3)
--small	Compressed fingerprint

### [General options]

-name	(prot) Name of molecule in out file	
-logf	Name of log file	
-D_Hb	N (3.5) Hbond length	(Angstroem)
-D_Hyd	N (4.5) Hydrophobic length	(Angstroem)
-D_Io	N (4.0) Ionic length	(Angstroem)
-D_Me	N (2.8) Metal/Acceptor length	(Angstroem)
-D_Ar	N (4.0) Aromatic interaction length	(Angstroem)
-D_Pic	N (4.0) Pi cation interaction length	(rad)
-a_H	N (Pi) HBond angle	(rad)
-at_H	N (Pi/3) HBond tolerance angle	(rad)
-a_ArFF	N (Pi) Aromatic Face to Face interaction angle	(rad)
-at_ArFF	N (Pi/6) Aromatic Face to Face tolerance angle	(rad)
-a_ArEF	N (Pi/2) Aromatic Edge to Face interaction angle	(rad)
-at_ArEF	N (Pi/3) Aromatic Edge to Face tolerance angle	(rad)
-a_Pic	N (Pi) Pi cation interaction angle	(rad)
-at_Pic	N (Pi/6) Pi cation tolerance angle	(rad)
--noMerge	Don't merge hydrophobic interactions	
--newH	Less permissive definition of hydrophobic interaction	

## GRIM - GRaph Interaction Matching :

<i>grim refProt refLig CompProt CompLig</i>	(1)
<i>grim refInts complInts</i>	(2)
<i>grim refProt refFile dockFile</i>	(3)

### [Note]

- (1) use --multim2 to use multimol2
- (3) refFile & dockFile can be multimol2 files

### [General options]

NOTE : INTERACTION GENERATION General Options also accessible

-rn	N (Ref) Reference name
-cn	N (Comp) Comparison name
--values	Only output score and not alignment
-sim	N (0) Boolean telling whether the pair is similar or not
-outInt	(MERG) Output only one kind of interaction positionning
MERG	All aligned interactions are outputed
LIG	InterLIG positionning
CENT	Centered positionning
PROT	InterPROT positionning

NOTE : outInt useless when used with --values

- match N (MERG) Align only with a specific position
- MERG Align with ALL interaction points
- LIG Align only with ligand interaction points
- PROT Align only with protein interaction points
- CENT Align only with centered interaction points
- score N (FCT) Scoring method function
- STD Scored by decreasing SumCI and increasing RMSD
- FCT Scored with scoring function
- newH Less permissive definition of hydrophobe interaction

[Alignment options]

- max N (1) Maximal number of outputed cliques.
- size N (3) Minimal size of a clique.
- all\_cliques Detect all cliques and not only maximal one
- dsame N (0.5) Maximal difference distance between to similar point(LIG-LIG, ...)
- dclose N (0.75) Maximal difference distance between to close point (LIG-CENT or PROT-CEN)
- dfar N (1) Maximal difference distance between to far point (LIG-CENT)

## VOLSITE - Cavity detection in a mol2 file

**vol site prot lig** (1)

**vol site prot** (2)

[General options]

- step N (1.5) Edge length of each box (Angstroem)
- boxS N (20) Edge length of the main box (Angstroem)
- b N (55) Minimal threshold for buriedness
- name N PDB Name for output cavity names
- n N (5) Minimal neighbours for buried cavity boxes
- nPTS N (70) Minimal number of cubes to consider it a cavity
- NPTS N (400) Maximal number of cubes to consider it a cavity
- dna Consider DNA as part of the protein
- cofactor Consider cofactor as part of the protein
- solvent Consider solvent as part of the protein
- hydrogen Consider hydrogens
- desc Write a descriptor file name descriptor.txt
- svm Build a svm property file
- drog N Observed druggability
- pharm Generate a pharmacophore (.chm) from cavity
- outExclu Output exclusion sphere in pharmacophore file

## PDB Process

**pdbconv protein[.pdb].mol2 output\_dir pdb\_id**

- wMOL2 Use MOL2 File as Input. PDB Options are not available
- wUnDrug Output undruggable cavities
- noLig PDB with no Ligand

By default all the following options are included. All chains will be kept

[PDB Options]

- HARMSIZE Harmonize size line to 80 characters
- MSEMET Change MSE to MET
- CSECYS Change CSE to CYS
- MOVHET move HETATM to the end of file
- ALTATM select alternative atoms
- NUMATM renumerotate atoms
- UPDMAS update the MASTER line

if you use one of the option below, you MUST use also --TOMOL2 option or use --wMOL2 option

--RESTYP	apply Residue Class (cofactor/STD_AA/MOD_AA/Ligand ...)
--BONDSE	create bonds
--CLNUNW	clean unwanted residues
--MOL2TY	apply MOL2 types according to templates
--SPLITM	split molecule into protein/ligand/solvent
--SelChain N	List of chains to keep, separated by underscore
--SELWAT	select water molecules
--SELLIG	select ligand

*utils bsa protein ligand*

*utils frag protein ligand*

*sims ref comp* (1)  
*sims file* (2)  
*sims Reflnt Complnt* (3)

--wints	To use interactions instead of fingerprints (option 3 only)
--small	Use small fingerprint (option 3 only)
--binary	To add when the fingerprint is binary
-metric N (TC)	Select the metric
TC	Tanimoto metric
HM	Hamming distance
RT	Ref Tversky
FT	Fit Tversky
DI	Dice
SO	Soergel

**detectppi** *PDB\_name name\_file* (1)

- c C : The name of all chains that might be in interaction
- n (1) : Maximum number of treated interfaces (ranked by decreasing size)

May have some bugs on structures with more than two chains

Z5&amp;24x3/V



## IChem Tutorial

### 1. License check and generation

The **IChem license** command enables you to see details of your license

---

#### > **IChem license**

---

The terminal output should look like this ...

```
Licence key generator : YES
Molecular realignment : YES
BSA Calculation : YES
IFP generator : YES
Interaction detection : YES
Graph Interaction Matching : YES
VolSite : YES
PDB to MOL2 conversion : YES
Patching MOL2 conversion : YES
Utils : YES
Sims : YES
Scoring : YES
Detection, analysis of PPI : YES
Your licence expires on (Y/M/D): 2019/12/31
```

For each module, a YES/NO flag indicate whether you have rights to use the corresponding utility. The last line gives the expiration date of your current license.

If your license scheme allows you to generate license keys, the **IChem** *genkey* command permits you to define any possible license file

**genKey** *year month day allowed\_tools*

year : Expiring year (4 digits : 2014)  
 month: Expiring month (2 digits : 09)  
 day : Expiring day (2 digits : 23)

allowed\_tools: each value is separated by space and must have a value of either 1 or 0. Must follow this order

First : Licence Key generator  
 Second : Realigning molecules  
 Third : BSA Calculation  
 Fourth : IFP generator  
 Fifth : Interaction detection  
 Sixth : Graph Interaction Matching  
 Seventh: VolSite  
 Eighth : PDB to MOL2 conversion  
 Ninth : Patching MOL2 conversion  
 Tenth : Utils  
 Eleventh: sims  
 Twelfth: scoring  
 Thirteenth: Detection, analysis of PPI

---

> **IChem** *genKey 2016 1 28 0 0 1 0 0 0 0 1 0 1 1*

---

grants you access until Jan.28<sup>th</sup> 2016 to 5 modules: BSA calculation, Patching MOL2 conversion, sims and scoring.

## 2. Protein-Ligand Interaction Fingerprints (IFPs)

The **IChem IFP** command registers in a bit string the interactions between a protein (active site) and a ligand.

**IFP protein ligand**

**IFP protein ligand ligand\_ref**

### [General options]

-name N (LIG) Name of the fingerprint -Default: Name of the ligand  
 --polar Detect and output only polar interactions  
 --metal Detect and output only metal interactions  
 --extended Include within the fingerprint:  
     |--> Metal/Acceptor interaction  
     |--> Weak Hydrogen bonds  
     |--> PI-Cation interactions

### [testing options]

-D\_Hb N (3.5) Hbond length (Angstroem)  
 -D\_Hyd N (4.5) Hydrophobic length (Angstroem)  
 -D\_Io N (4.0) Ionic length (Angstroem)  
 -D\_Me N (2.8) Metal/Acceptor length (Angstroem)  
 -D\_Ar N (4.0) Aromatic interaction length (Angstroem)  
 -D\_Pic N (4.0) Pi cation interaction length (rad)  
 -a\_H N (Pi) HBond angle (rad)  
 -at\_H N (Pi/3) HBond tolerance angle (rad)  
 -a\_ArFF N (Pi) Aromatic Face to Face interaction angle (rad)  
 -at\_ArFF N (Pi/6) Aromatic Face to Face tolerance angle (rad)  
 -a\_ArEF N (Pi/2) Aromatic Edge to Face interaction angle (rad)  
 -at\_ArEF N (Pi/3) Aromatic Edge to Face tolerance angle (rad)  
 -a\_Pic N (Pi) Pi cation interaction angle (rad)  
 -at\_Pic N (Pi/6) Pi cation tolerance angle (rad)  
 --ligD Print all possible ligand interactions

Please note that ligand file can be a multi-mol2 file

For more information, see:



Marcou G, Rognan D.

J Chem Inf Model. 2007 Jan-Feb;47(1):195-207.

[Optimizing fragment and scaffold docking by use of molecular interaction fingerprints.](#)

In a first mode, interactions between an active site (mol2 file format) and one/several ligands (single or multimol2 file) are outputted in a table along with the IFP bit string. In a second mode, similarity of the IFP(s) to that of one (several) references is outputted in addition.

**Output file:** **ligand.ifp**

```
> IChem IFP site.mol2 ligand.mol2 >ligand.ifp
```

HBond_LIG	OD1	78	ASP 113-A	O17	1	CAU	0-XX	0	2.60795	134.667
HBond_PROT	ND2	421	ASN 312-A	O17	1	CAU	0-XX	1	2.76961	159.236
Hydrophobic	CZ	349	PHE 289-A	C16	2	CAU	0-XX	3	3.74597	/
Ionic_LIG	OD2	79	ASP 113-A	N19	4	CAU	0-XX	4	2.94042	/
Ionic_LIG	OD1	78	ASP 113-A	N19	4	CAU	0-XX	5	3.60381	/
Hydrophobic	CZ3	46	TRP 109-A	C21	6	CAU	0-XX	6	4.43435	/
Hydrophobic	CG2	64	THR 110-A	C21	6	CAU	0-XX	7	4.42961	/
Hydrophobic	CB	163	PHE 193-A	C21	6	CAU	0-XX	8	4.07199	/
Hydrophobic	CH2	47	TRP 109-A	C22	7	CAU	0-XX	9	3.72321	/
Hydrophobic	CZ	438	TYR 316-A	C22	7	CAU	0-XX	10	4.46257	/
Hydrophobic	CG2	125	VAL 117-A	C15	8	CAU	0-XX	11	4.0788	/
Hydrophobic	CE2	348	PHE 289-A	C15	8	CAU	0-XX	13	4.0318	/
HBond_LIG	OG	248	SER 203-A	N7	16	CAU	0-XX	23	3.31671	127.808
Hydrophobic	CG2	185	THR 195-A	C1	20	CAU	0-XX	24	4.32061	/
Hydrophobic	CZ3	327	TRP 286-A	C16	2	CAU	0-XX	25	4.05558	/
Hydrophobic	CG1	124	VAL 117-A	C12	11	CAU	0-XX	28	3.98025	/
Hydrophobic	CG1	89	VAL 114-A	C11	12	CAU	0-XX	29	3.83646	/
Hydrophobic	CB	288	SER 207-A	C10	13	CAU	0-XX	30	3.68825	/
Hydrophobic	CG2	90	VAL 114-A	C13	10	CAU	0-XX	31	3.81147	/

7<sup>th</sup> position: ionic (ligand charged +)

Page 15

**Input directory:** \$ICHEM\_DIR/test/IFP

**Input files:** site.mol2, docked.mol2, ligand.mol2

**Output file:** docked Tc.ifp

```
> IChem IFP site.mol2 docked.mol2 ligand.mol2 >docked Tc.ifp
```

The following information is contained in the output file (docked\_Tc.ifp) or at the terminal if no file redirection (*... >...*) is given:

REF |A M82|A V86|A W109|A T110|A D113|A V114|A L115|A V117|A T118|A C191|A F193|A T195|A  
Y199|A A200|A I201|A S203|A S204|A I205|A S207|A F208|A W286|A F289|A F290|A N293|A Y308|A  
N312|A Y316

```
REF
00000000000000100000010000000001011000000000000100000000000000000000000100000010000000000000000
000000000000000010000000000000001000000000000100000010000000000000000000000000000000000000000000000
0
```

DOCKED |A M82|A V86|A W109|A T110|A D113|A V114|A L115|A V117|A T118|A C191|A F193|A T195|A  
Y199|A A200|A I201|A S203|A S204|A I205|A S207|A F208|A W286|A F289|A F290|A N293|A Y308|A  
N312|A Y316

DOCKED

000000000000000100000010000001100000000000001000000000000000000100000010000000000000000

000000000000000000000000000000010000000000000100000010000000000000100000010000000001000100000

0

REF	REF	1
REF	DOCKED	0.705882

In the first section, IFPs are just given for the reference (REF) and the docked ligand (DOCKED)

In the second section, the Tc coefficients are given for both the REF and the DOCKED ligand

### -Non standard modes

Polar mode: only 5 bits/residues are outputted (positions 4-7 of the standard IFP + metal coordination)

> IChem --polar IFP arguments

Extended mode: 2 additional bits/residues with respect to the normal mode (position 8: pi-cation interactions; position 9: metal coordination)

> IChem --extended IFP arguments

### -Interaction detection rules

-D_Hb	N (3.5)	Hbond length	(Angstroem)
-D_Hyd	N (4.5)	Hydrophobic length	(Angstroem)
-D_lo	N (4.0)	Ionic length	(Angstroem)
-D_Me	N (2.8)	Metal/Acceptor length	(Angstroem)
-D_Ar	N (4.0)	Aromatic interaction length	(Angstroem)
-D_Pic	N (4.0)	Pi cation interaction length	(Angstroem)
-a_H	N (Pi)	HBond angle	(rad)
-at_H	N (Pi/3)	HBond tolerance angle	(rad)
-a_ArFF	N (Pi)	Aromatic Face to Face interaction angle	(rad)
-at_ArFF	N (Pi/6)	Aromatic Face to Face tolerance angle	(rad)



-a\_ArEF      N (Pi/2) Aromatic Edge to Face interaction angle      (rad)  
 -at\_ArEF    N (Pi/3) Aromatic Edge to Face tolerance angle      (rad)  
 -a\_Pic       N (Pi)    Pi cation interaction angle                      (rad)  
 -at\_Pic       N (Pi/6) Pi cation tolerance angle                      (rad)  
 --ligD       Print all possible ligand interactions

Interactions are detected based on default topological rules that can be modified whenever necessary:

> **IChem** *-D\_Hb DHB -D\_Hyd DHYD -D\_Io DIO -D\_Me DME -D\_Ar DAR -D\_PIC DIC -a\_H AH -at\_H ATH -a\_ArFF AARFF -at\_ArFF ATARFF -a\_ArEF -AAREF -at\_AREF ATAREF -a\_Pic APIC -at\_PIC ATPIC* **IFP arguments**

Default values are described in the following table:

Parameter	Description	Default Value
DHB	Donor-acceptor distance threshold for H-bond interaction	3.5 Å
DHYD	Donor-acceptor distance threshold for hydrophobic interaction	4.5 Å
DIO	Anion-cation distance threshold for ionic interaction	4.0 Å
DME	Metal-acceptor distance threshold for metal chelation	2.8 Å
DAR	Distance threshold between two aromatic ring centers	4.0 Å
DPIC	Distance threshold for pi cation interaction	4.0 Å
AH	Angle threshold between donor-H-acceptor for H-bond	$\pi$
ATH	Tolerance angle for H-bond	$\pi/3$
AARFF	Angle threshold between the plane of two aromatic cycles for face to face aromatic interactions	$\pi$
ATARFF	Tolerance angle for face to face aromatic interactions	$\pi/6$
AAREF	Angle threshold between the plane of two aromatic cycles for edge to face aromatic interactions	$\pi/2$
ATAREF	Tolerance angle for edge to face aromatic interactions	$\pi/6$
APIC	Angle threshold for pic cation interactions	$\pi$
ATPIC	Tolerance angle for pi cation interactions	$\pi/6$
nomerge	After detecting all hydrophobic contacts, a filtering process is engaged to avoid too many hydrophobic pseudoatoms. Using this flag avoids the filtering process	no

The interaction rules options can be used for the following IChem tools: **IFP**, **ints**, **grim**



**Some notes:**

- Beware: the mol2 file should comply with the standard mol2 file format from TRIPOS (<http://www.tripos.com/data/support/mol2.pdf>)
- If you do not have SYBYL, you can generate correct mol2 files with UCSF chimera (<http://www.cgl.ucsf.edu/chimera>)
- IFPs are active-site dependent (nbits/residue) and cannot be compared across active sites of different compositions (e.g. different number of residues).
- Please rather use the active site and not the protein as input, to restrict the size of the IFP and prevent generating bit strings with mostly "0" values.
- Compare your output with that given in \$ICHEM\_DIR/test/IFP/output

### 3. Interaction Fingerprint Triplet (TIFPs)

TIFPs are coordinate-frame invariant interaction fingerprints presenting the advantage to compare ligand binding to completely different active sites, whatever their composition.


For more information, see:



Desaphy J, Raimbaud E, Ducrot P, Rognan D.  
J Chem Inf Model. 2013 Mar 25;53(3):623-37

[Encoding protein-ligand interaction patterns in fingerprints and graphs.](#)

Interactions are encoded by pseudoatoms in a MOL2 file format with the following properties: interaction type, atomic coordinates (**Figure 1**).



**Figure 1.** Encoding protein-ligand interaction by pseudoatoms. Each interaction is detected on the fly from standard topological rules and is represented by a pseudoatom with 3 possible atomic coordinates: the interacting ligand atom (LIG mode, violet balls), the interacting protein atom (PROT mode, orange balls), the barycenter of interacting protein and ligand atoms (CENT mode), the three aboved-described atoms (MERG mode).

### ints prot lig out

```
ints prot lig outfile
  -type          (CENT) Alter positionning output
                    Multiple values allowed, separated by space
                    PROT      InterPROT positionning
                    LIG       InterLIG positionning
                    CENT      Centered positionning
                    MERG      Merged all 3 above

  -fgps          (STD)  Fingerprint format
                    STD       Standard (1 0 21 0 0 3)
                    SVM       SVM format (1:1 3:21 6:3)
                    CMP       Compressed (1 [1 21 [2 3)
  --small        Compressed fingerprint

[General options]
  -name          (prot) Name of molecule in out file
  -logf          Name of log file
  -D_Hb          N (3.5) Hbond length (Angstroem)
  -D_Hyd         N (4.5) Hydrophobic length (Angstroem)
  -D_Io          N (4.0) Ionic length (Angstroem)
  -D_Me          N (2.8) Metal/Acceptor length (Angstroem)
  -D_Ar          N (4.0) Aromatic interaction length (Angstroem)
  -D_Pic         N (4.0) Pi cation interaction length (Angstroem)
  -a_H           N (Pi) HBond angle (rad)
  -at_H          N (Pi/3) HBond tolerance angle (rad)
  -a_ArFF        N (Pi) Aromatic Face to Face interaction angle (rad)
  -at_ArFF       N (Pi/6) Aromatic Face to Face tolerance angle (rad)
  -a_ArEF        N (Pi/2) Aromatic Edge to Face interaction angle (rad)
  -at_ArEF       N (Pi/3) Aromatic Edge to Face tolerance angle (rad)
  -a_Pic         N (Pi) Pi cation interaction angle (rad)
  -at_Pic        N (Pi/6) Pi cation tolerance angle (rad)
  --noMerge      Don't merge hydrophobic interactions
  --newH         Less permissive definition of hydrophobic interaction
```

### -Listing interactions and outputting interaction pseudoatoms

Input directory: \$ICHEM\_DIR/test/TIFP

Input files: site.mol2, ligand.mol2

Output files: 2rh1\_ints.txt, 2rh1\_site\_INTS\_C.mol2

---

```
> IChem -logf 2rh1_ints.txt -type CENT ints site.mol2 ligand.mol2
```

---

The 2rh1\_ints.txt file contains the same information as the table outputted by the IChem IFP command (recall pages 13, 14)

The 2rh1\_site\_INTS\_X.mol2 file describes the properties and atomic coordinates of the interaction pseudoatoms. Depending on the -type option (CENT, PROT, LIG, MERG), pseudoatoms are mapped onto ligand-interacting atoms (-type LIG), protein-interacting atoms (-type PROT), barycenter of ligand and protein-interacting atoms (-type CENT, default option), or on all above-described atoms (-type MERG)

Filename	-type mode
Header_INTS_C.mol2	CENT mode
Header_INTS_L.mol2	LIG mode
Header_INTS_P.mol2	PROT mode
Header_INTS_M.mol2	MERG mode

### Atomic description of the interaction pseudoatoms:

Atom Name	Described interaction
CA	hydrophobic
CZ	aromatic
O	Hydrogen-bond (interacting protein atom is acceptor)
OG	Hydrogen-bond (interacting protein atom is both acceptor and donor)
N	Hydrogen-bond (interacting protein atom is donor)
OD1	Ionic (interacting protein atom is negatively charged)
NZ	Ionic (Interacting protein atom is positively charged)
ZN	Metal coordination

### -Generating interaction fingerprints

With respect to *IChem IFP*, *IChem ints* outputs generic and binding site-independent interaction fingerprints consisting in all possible combinations of triplets of interaction pseudoatoms

Input directory: *\$ICHEM\_DIR/test/TIFP*

Input files: *site.mol2*, *ligand.mol2*

Output files: *2rh1.fgp*

---

> *IChem -fgps STD ints site.mol2 ligand.mol2 2rh1\_full.fgp*

---

The fgp output file can be appended after every novel *IChem ints* command if the same outputfile name is given. By default, a standard fingerprint format is outputted, but the *-fgps* option enables you to control the fingerprint format:

Format mode	<i>-fgps option</i>	Example	Description
Standard	STD	0 1 0 0 132 0 1 0 0 0	Each value is separated by a space Inefficient for long and sparse fingerprints (many null values)
SVM	SVM	2:1: 5:132 7:1	Only non-null values are outputted and defined by <i>POSITION:VALUE</i> separated by a space Very efficient for long sparse fingerprints Not efficient for long and dense fingerprints (few null values)
Compressed	CMP	[1 1 [2 132 [1 1 [3	Non-null values are explicitly outputted. Null values are encoded by a "[" sign followed by the number of consecutive null values Efficient for long and sparse fingerprints

By default, the standard fingerprint is made of 12510 integers, most of which encoding triplets of pseudoatoms describing hydrophobic contacts. The full fingerprint can be pruned to remove rare triplets (based on the analysis of ca. 10 000 protein-ligand PDB Structures) and lead to a simpler fingerprint of 210 integers (*--small* option)

---

> *IChem --small -fgps STD ints site.mol2 ligand.mol2 2rh1\_small.fgp*

---

## 4. Graph Matching (GRIM)

GRIM is a tool to match protein-ligand complexes using a graph matching algorithm focusing on protein-ligand interaction pseudoatoms. It thereby enables an interaction-based alignment of different protein-ligand complexes that can be quantified by an empirical scoring function (GrScore) and used to post-process docking poses by similarity to known protein-ligand interaction patterns.

In a recent international contest (D3R Docking Challenge 2015), **Grim** was ranked **2<sup>nd</sup> out of 44** scoring functions to predict the binding mode of 36 inhibitors prior to the release of protein-bound X-ray coordinates ([Gathiaka et al. J Comput.-Aided Mol. Des, 2016](#))

For more information, see:



Desaphy J, Raimbaud E, Ducrot P, Rognan D., J Chem Inf Model. 2013; 53: 623-637:  
[Encoding protein-ligand interaction patterns in fingerprints and graphs.](#)

Slyenko I, Da Silva F, Bret G, Rognan D. J Comput Aided Mol Des. 2016, 30: 669-683.  
[Docking pose selection by interaction pattern graph similarity: application to the D3R grand challenge 2015.](#)

<b>grim refProt refLig CompProt CompLig</b>	<b>(1)</b>
<b>grim refInts complnts</b>	<b>(2)</b>
<b>grim refProt refFile dockFile</b>	<b>(3)</b>

[Note]

- (1) use --multim2 to use multimol2 ligand files
- (3) refFile & dockFile can be multimol2 files

[General options]

NOTE : INTERACTION GENERATION General Options also accessible

```
-rn      N (Ref)   Reference name
-cn      N (Comp)  Comparison name
--values          Only output score and not alignment
-sim      N (0)    Boolean telling whether the pair is similar or not
-outInt     (MERG) Output only one kind of interaction positioning
    MERG          All aligned interactions are outputed
    LIG           InterLIG positioning
    CENT          Centered positioning
    PROT          InterPROT positioning
    NOTE : outInt useless when used with --values
-match  N (MERG)   Align only with a specific position
    MERG          Align with ALL interaction points
    LIG           Align only with ligand interaction points
    PROT          Align only with protein interaction points
    CENT          Align only with centered interaction points
-score  N (FCT)    Scoring method function
    STD          Scored by decreasing SumCl and increasing RMSD
    FCT          Scored with scoring function
--newH          Less permissive definition of hydrophobe interaction
```

[Alignment options]

```
-max      N (1)    Maximal number of outputed cliques.
-size     N (3)    Minimal size of a clique.
--all_cliques  Detect all cliques and not only maximal one
-dsame    N (0.5)  Maximal difference distance between to similar point (LIG-
LIG, ...)
-dclose   N (0.75) Maximal difference distance between to close point (LIG-
CENT or PROT-CENT)
-dfar     N (1)    Maximal difference distance between to far point (LIG-CENT)
-1st mode: Graph matching from structures
```

**Input directory:** \$ICHEM\_DIR/test/GRIM

**Input files:** 2rh1\_prot.mol2, 2rh1\_lig.mol2, 4amj\_prot.mol2, 4amj\_lig.mol2

**Output files:** Complnts.mol2, Grifp\_res.csv, GRIM\_ints.mol2, GRIM\_lig.mol2, GRIM.log,  
GRIM\_prot.mol2, RefInts.mol2

---

```
>IChem -sim 1 -rn 2rh1 -cn 4amj grim 2rh1_prot.mol2 2rh1_lig.mol2 4amj_prot.mol2
4amj_lig.mol2
```

---

**-2<sup>nd</sup> mode: Graph matching from interaction pseudoatoms**

**Input directory:** \$ICHEM\_DIR/test/GRIM

**Input files:** 2rh1\_INTS\_M.mol2, 4amj\_INTS\_M.mol2

**Output files:** Complnts.mol2, Grifp\_res.csv, GRIM\_ints.mol2, GRIM\_lig.mol2, GRIM.log, GRIM\_prot.mol2, Reflnts.mol2

---

```
> IChem -sim 1 -rn 2rh1 -cn 4amj grim 2rh1_INTS_M.mol2 4amj_INTS_M.mol2
```

---

Whatever the mode, 7 files are outputted:

**Complnts.mol2:** interaction pseudoatoms (merged mode) of the complex to fit

**Grifp\_res.csv:** summary of the GRIM alignment

**GRIM\_ints.mol2:** aligned interaction pseudoatoms of the complex to fit

**GRIM\_lig.mol2:** coordinates of the fitted ligand, aligned to the reference

**GRIM.log:** output of the GRIM alignment

**GRIM\_prot.mol2:** coordinates of the fitted protein, aligned to the reference

**Reflnts.mol2:** interaction pseudoatoms (merged mode) of the reference complex

The **Grifp\_res.csv** file (appended at every novel IChem grim command) is a table that looks like this:

NCli	Ref	Comp	Simil	LIG	CENTER	PROT	SumCl	RMSD	RLig	RCent	RProt	CLig	CCent	CProt	GrSc	RMSDAI	NPol
1	2rh1	4amj	1	9	6	7	0.0685	0.2763	13	19	19	15	21	21	0.75450	0.51	
589	10																

NCli: Id of the clique (starts at 0)

Ref: Name of the reference (protein name if **-rn** is not used)

Comp: Name of the comparison (protein name if **-cn** is not used)

Simil: Similarity flag (empty is **-sim** is not used)

LIG: Number of matched LIG pseudoatoms

CENTER: Number of mgrimatched CENTERED pseudoatoms

PROT: Number of matched PROT pseudoatoms

SumCl: SumCl score (clique-based)

RMSD: root-mean square deviation (in Å) of the clique

RLig: Number of reference LIG pseudoatoms

RCent: Number of reference CENTERED pseudoatoms

RProt: Number of reference PROT pseudoatoms

CLig: Number of comparison LIG pseudoatoms

CCent: Number of comparison CENTERED pseudoatoms

CProt: Number of comparison PROT pseudoatoms

GrSC: Grim score (empirical score)

RMSDAI: root-mean square deviation (in Å) of the clique according to GrScore

NPol: Number of matched polar interaction pseudoatoms

**-3<sup>rd</sup> mode: Post-processing docking results by interaction graph matching**

**Input directory:** \$ICHEM\_DIR/test/GRIM/screen

**Input files:** docked.mol2, ligand.mol2, site.mol2

**Output files:** Grimscreen\_res.tsv, GrScreen.csv

> IChem -sim 1 grim site.mol2 ligand.mol2 docked.mol2

The file docked.txt summarized the docking results obtained with Surflex-Dock: pose nr, docking score (pkd), rmsd to the true X-ray pose

Pose	pkd	rmsd
ligand_000	8.63	1.14
ligand_001	8.35	1.08
ligand_002	8.12	1.00
ligand_003	8.07	1.14
ligand_004	7.54	1.19
ligand_005	7.10	1.28
ligand_006	7.04	1.41
ligand_007	6.90	1.08
ligand_008	6.72	1.17
ligand_009	6.67	1.34
ligand_010	6.46	1.23
ligand_011	6.34	1.47
ligand_012	6.31	2.02
ligand_013	6.30	3.67
ligand_014	6.28	1.50
ligand_015	6.28	1.27
ligand_016	6.20	2.09
ligand_017	6.16	1.18
ligand_018	6.07	3.07
ligand_019	6.04	3.54

All docking poses (docked.mol2) are matched to the X-ray pose (ligand.mol2) for rescoring based on the interaction pattern graph similarity score (GrScore). The full output is stored in the Grimscreen\_res.tsv output file.

NCLi	Ref	Comp	Simil	LIG	CENTER	PROT	SumCl	RMSD	RLig	RCent	RProt	CLig	CCent	CProt	GrSc	RMSDAL	NPol
1	2rhl_site	REF	1	9	13	13	0.1054	0.2403	13	19	19	13	19	18	0.8779	0.0049	11
1	2rhl_site	REF	1	10	11	11	0.0924	0.1405	13	19	19	13	18	17	0.8611	0.0023	9
1	2rhl_site	REF	1	11	14	16	0.1151	0.2133	13	19	19	14	21	21	0.9589	0.0026	14
1	2rhl_site	REF	1	10	12	13	0.1064	0.2206	13	19	19	14	20	20	0.8907	0.0024	13
1	2rhl_site	REF	1	8	9	5	0.0703	0.2583	13	19	19	13	17	16	0.7242	0.0075	7
1	2rhl_site	REF	1	9	6	4	0.0538	0.2686	13	19	19	12	18	16	0.7091	0.0151	5
1	2rhl_site	REF	1	9	5	5	0.0565	0.2744	13	19	19	13	20	17	0.7196	0.0186	6
1	2rhl_site	REF	1	7	12	16	0.1199	0.2240	13	19	19	12	17	16	0.8909	0.0039	10
1	2rhl_site	REF	1	11	11	13	0.1047	0.3423	13	19	19	14	23	22	0.8928	0.0030	14
1	2rhl_site	REF	1	8	4	7	0.0573	0.4584	13	19	19	13	19	17	0.7161	0.0072	6
1	2rhl_site	REF	1	5	6	6	0.0527	0.3720	13	19	19	14	21	18	0.6696	0.0075	7
1	2rhl_site	REF	1	4	2	6	0.0313	0.4146	13	19	19	12	17	16	0.6313	0.0150	2
1	2rhl_site	REF	1	3	3	3	0.0191	0.5820	13	19	19	10	16	15	0.5626	0.2846	0
1	2rhl_site	REF	1	4	3	6	0.0303	0.5961	13	19	19	13	19	17	0.6219	0.0020	2
1	2rhl_site	REF	1	4	3	6	0.0396	0.4205	13	19	19	12	17	15	0.6365	0.0020	3
1	2rhl_site	REF	1	7	5	9	0.0665	0.3336	13	19	19	13	16	15	0.7444	0.0031	5
1	2rhl_site	REF	1	5	8	6	0.0427	0.4845	13	19	19	12	19	18	0.6672	0.0300	2
1	2rhl_site	REF	1	6	11	15	0.0853	0.2214	13	19	19	14	19	17	0.8506	0.0033	7
1	2rhl_site	REF	1	4	1	2	0.0128	0.3104	13	19	19	12	19	15	0.5738	0.1752	0
1	2rhl_site	REF	1	4	7	7	0.0341	0.4100	13	19	19	13	20	19	0.6660	0.0043	1

The pose with the highest Grimscore is not the top-ranked one (ligand\_000) but pose 002 (highlighted in yellow). Interestingly, it is the pose with the smallest rmsd to the true X-ray pose.



## - GRIM-based alignment of protein-ligand complexes

*realign rigidM mobilM applied1 applied2*

||-> rigidM : reference molecule to apply alignment to  
 ||-> mobilM : comparison molecule to apply alignment from  
 ||-> applied: molecule to apply rotation/translation to

[General options]

-gmatch N (NAME) Use graph matching to align  
 NAME Atom Name matching  
 ATMN Atomic Name matching  
 MOL2 MOL2 Type matching  
 CALP CAlpha Atom matching (protein only)  
 --wMob Also outputs the aligned mobilM  
 -rule R

By default, the program will perform an atom by atom match, without taking care of what kind of atom it match. If you want to perform a match by regarding only some atoms, this index\_string is here to do so

ex : -i '2-3|1-6|23-160' Will match the second atom from the reference with the third from the comparison, the first with the sixth ...

Input directory: **\$ICHEM\_DIR/test/REALIGN**

Input files: **4amj\_INTS\_M.mol2, GRIM\_ints.mol2, 4amj\_prot.mol2 4amj\_lig.mol2**

Output files: **rot\_4amj\_INTS\_M.mol2, rot\_4amj\_lig.mol2, rot\_4amj\_prot.mol2**

---

```
> IChem --wMob realign GRIM_ints.mol2 4amj_INTS_M.mol2 4amj_lig.mol2
4amj_prot.mol2
```

---

**IChem *realign*** will apply to the native ligand and/or protein input files (4amj\_lig.mol2, 4amj\_prot.mol2) a set of rotation/translations (deduced from the previous step) to align them to the reference GRIM input files (2rh1\_lig.mol2, 2rh1\_prot.mol2). The move mimics the transformation of interaction pseudoatoms previously saved: **4amj\_INTS\_M.mol2 → GRIM\_ints.mol2**

The aligned coordinates of the complex to fit are in the rot\_xxx.mol2 files.

## 5. Cavity detection and druggability prediction (VolSite)

**IChem VolSite** is a structure-based tool to automatically detect cavities at the surface of a target protein, and predict their ligandability (structural druggability)

For more information, see:



Desaphy J, Azdimousa K, Kellenberger E, Rognan D.

J Chem Inf Model. 2012 Aug 27;52(8):2287-99.

[Comparison and druggability prediction of protein-ligand binding sites from pharmacophore-annotated cavity shapes.](#)

A new feature, introduced from the current version is the generation of cavity-based pharmacophores, as described in

Tran-Nguyen VK, Da Silva F, Bret G, Rognan D.

J Chem Inf Model. 2019, 59, 1, 573-585

[All in One: Cavity Detection, Druggability Estimate, Cavity-Based Pharmacophore Perception and Virtual Screening.](#)

**IChem Volsite** can be run in two modes depending on whether coordinates of a bound ligand are given (ligand-restricted mode) or not (unrestricted mode).

**volsite prot lig**

**(1)**

**volsite prot**

**(2)**

```
[General options]
-step N (1.5)      Edge length of each box (Angstroem)
-boxS N (20)       Edge length of the main box (Angstroem)
-b N (55)          Minimal threshold for buriedness
-name N           PDB Name for output cavity names
-n N (5)           Minimal neighbours for buried cavity boxes
-nPTS N (70)       Minimal number of cubes to consider it a cavity
-NPTS N (400)      Maximal number of cubes to consider it a cavity
--dna             Consider DNA as part of the protein
--cofactor        Consider cofactor as part of the protein
--solvent         Consider solvent as part of the protein
--hydrogen        Consider hydrogens
--desc            Write a descriptor file name descriptor.txt
--svm             Build a svm property file
--drog N          Observed druggability
--pharm           Generate a pharmacophore (.chm) from cavity
--outExclu        Output exclusion sphere in pharmacophore file
```

### -Detection of all possible cavities (unrestricted mode)

Input directory: **\$ICHEM\_DIR/test/Volsite/unrestricted**

Input files: **protein.mol2**

Output files: **CAVITY\_Nx\_ALL.mol2 (x=1-9), VolSite\_stat.csv**

---

**> IChem Volsite protein.mol2**

---

VolSite produces a summary (**VolSite\_stat.csv**) and a mol2 file (**CAVITY\_Nx\_ALL.mol2; x=1-6**) for each detected cavity. For the test example (**protein.mol2**), 6 cavities are detected at the surface of the protein and numbered (1 to 6) from the largest to the smallest.

The **VolSite\_stat.csv** summary recapitulates the properties of all detected cavities as follows:

Name	NCav	Size	Buriedness	NPts	Volume	CA	CZ	O	OG	OD1	N	NZ	DU	Drugg
/	1	ALL	84.9875	321	1083.38	87	76	11	27	8	61	34	17	1.12574
/	2	ALL	78.1878	149	502.875	63	15	4	15	8	12	13	19	0.71907
/	3	ALL	70.0841	133	448.875	30	9	6	5	0	32	28	23	-0.0102194
/	4	ALL	71.5261	80	270	36	11	3	5	2	7	0	16	0.543419
/	5	ALL	71.0452	78	263.25	53	9	0	0	0	2	0	14	1.05675
/	6	ALL	71.3473	75	253.125	18	5	5	0	18	1	18	10	-1.98361

1<sup>st</sup> column: name of the protein (by default, mol2 file protein header)

2<sup>nd</sup> column: Cavity number

3<sup>rd</sup> column: Size (ALL: no distance cut-off to the cavity center is applied to define cavity points)

4<sup>th</sup> column: Average buriedness of the cavity, in %

6<sup>th</sup> column: Cavity volume in Å<sup>3</sup>

7<sup>th</sup> column: number of hydrophobic cavity points (CA)

8<sup>th</sup> column: number of aromatic cavity points (CZ)

9<sup>th</sup> column: number of hydrogen-bond acceptor cavity points (O)

10<sup>th</sup> column: number of hydrogen-bond acceptor/donor cavity points (OG)

11<sup>th</sup> column: number of negatively ionizable cavity points (OD1)

12<sup>th</sup> column: number of hydrogen bond donor cavity points (N)

13<sup>th</sup> column: number of positively ionized cavity points (NZ)

14<sup>th</sup> column: number of dummy cavity points (DU): no protein atom < 4.5 Å

15<sup>th</sup> column: estimated druggability (Drugg): druggable if Drugg > 0; undruggable if Drugg < 0

#### -Detection of a cavity around a particular ligand (ligand-restricted mode)

Input directory: \$ICHEM\_DIR/test/Volsite/ligand

Input files: **protein.mol2**, **ligand.mol2**

Output files: **CAVITY\_Nx\_y.mol2** (x = 1-3; y = 4, 6, 8, 12, ALL), **VolSite\_stat.csv**

---

> IChem **volsite** **protein.mol2** **ligand.mol2**

---

The terminal output of the command indicates that a druggable cavity has been found (Druggability: 1.13).

VolSite produces a summary (**VolSite\_stat.csv**) and a mol2 file (**CAVITY\_Nx\_y.mol2**; x = 1; y = 4, 6, 8, 12, ALL) for each detected cavity. For the test example (**protein.mol2**, **ligand.mol2**), 2 cavities are detected at the surface of the protein and numbered (x= 1 to 2) from the largest to the smallest. y indicates the maximal distance (in Å) between cavity points and any input ligand heavy atom. For each cavity, 5 truncation modes are then applied to define binding sites of increasing size (4, 6, 8, 10, 12 Å) around the input ligand. If y = ALL, no truncation is applied.

The **VolSite\_stat.csv** summary recapitulates the properties of all detected cavities as follows:

Name	NCav	Size	NPts	Volume	CA	CZ	O	OG	OD1	N	NZ	DU	Recovery	Drugg
CAU	1	4	115	388.125	31	37	3	12	0	22	6	4	53.913	/
CAU	1	6	160	540	44	44	6	20	0	30	10	6	38.75	/
CAU	1	8	205	691.875	58	56	6	22	1	37	15	10	30.2439	/
CAU	1	12	303	1022.62	83	74	9	26	8	56	30	17	20.462	/
CAU	1	ALL	321	1083.38	87	76	11	27	8	61	34	17	19.3146	1.12574

2

1<sup>st</sup> column: name of the ligand (by default, mol2 file ligand header)

2<sup>nd</sup> column: Cavity number

3<sup>rd</sup> column: Size (truncation mode)

4<sup>th</sup> column: Number of cavity points

5<sup>th</sup> column: Cavity volume in Å<sup>3</sup>  
 6<sup>th</sup> column: number of hydrophobic cavity points (CA)  
 7<sup>th</sup> column: number of aromatic cavity points (CZ)  
 8<sup>th</sup> column: number of hydrogen-bond acceptor cavity points (O)  
 9<sup>th</sup> column: number of hydrogen-bond acceptor/donor cavity points (OG)  
 10<sup>th</sup> column: number of negatively ionizable cavity points (OD1)  
 11<sup>th</sup> column: number of hydrogen bond donor cavity points (N)  
 12<sup>th</sup> column: number of positively ionized cavity points (NZ)  
 13<sup>th</sup> column: number of dummy cavity points (DU): no protein atom < 4.5 Å  
 14<sup>th</sup> column: % of the binding site enclosing the ligand  
 15<sup>th</sup> column: estimated druggability (Drugg): druggable if Drugg > 0; undruggable if Drugg < 0

### - Outputting cavity properties

The **IChem --desc** command outputs cavity properties used by our SVM druggability predictor.

**Input directory:** \$ICHEM\_DIR/test/Volsite/unrestricted

**Input files:** protein.mol2

**Output files:** CAVITY\_Nx\_ALL.mol2 (x=1-6), VolSite\_stat.csv, 2rh1\_protein\_descriptor.txt

---

```
> IChem --desc volsite protein.mol2
```

---

The 2rh1\_protein\_descriptor.txt file outputs for every cavity a vector of 89 reals as follows:

#	Descriptor
1	Volume
2	% of aromatic cavity points (CZ)
3	% of hydrophobic points (CA)
4	% of h-bond acceptor points (O)
5	% of negatively ionizable points (OD1)
6	% of h-bond acceptor & donor points (OG)
7	% of h-bond donor points (N)
8	% of positively ionizable points (NZ)
9	% of dummy points (DU)
10	Percent of CZ points with a projection value below 40
11	Percent of CZ points with a projection value between 40 and 50
12	Percent of CZ points with a projection value between 50 and 60
13	Percent of CZ points with a projection value between 60 and 70
14	Percent of CZ points with a projection value between 70 and 80
15	Percent of CZ points with a projection value between 80 and 90
16	Percent of CZ points with a projection value between 90 and 100
17	Percent of CZ points with a projection value between 100 and 110
18	Percent of CZ points with a projection value between 110 and 120
19	Percent of CZ points with a projection value of 120
20	Percent of CA points with a projection value below 40
21	Percent of CA points with a projection value between 40 and 50
22	Percent of CA points with a projection value between 50 and 60
23	Percent of CA points with a projection value between 60 and 70
24	Percent of CA points with a projection value between 70 and 80
25	Percent of CA points with a projection value between 80 and 90

26	Percent of CA points with a projection value between 90 and 100
27	Percent of CA points with a projection value between 100 and 110
28	Percent of CA points with a projection value between 110 and 120
29	Percent of CA points with a projection value of 120
30	Percent of O points with a projection value below 40
31	Percent of O points with a projection value between 40 and 50
32	Percent of O points with a projection value between 50 and 60
33	Percent of O points with a projection value between 60 and 70
34	Percent of O points with a projection value between 70 and 80
35	Percent of O points with a projection value between 80 and 90
36	Percent of O points with a projection value between 90 and 100
37	Percent of O points with a projection value between 100 and 110
38	Percent of O points with a projection value between 110 and 120
39	Percent of O points with a projection value of 120
40	Percent of OD1 points with a projection value below 40
41	Percent of OD1 points with a projection value between 40 and 50
42	Percent of OD1 points with a projection value between 50 and 60
43	Percent of OD1 points with a projection value between 60 and 70
44	Percent of OD1 points with a projection value between 70 and 80
45	Percent of OD1 points with a projection value between 80 and 90
46	Percent of OD1 points with a projection value between 90 and 100
47	Percent of OD1 points with a projection value between 100 and 110
48	Percent of OD1 points with a projection value between 110 and 120
49	Percent of OD1 points with a projection value of 120
50	Percent of OG points with a projection value below 40
51	Percent of OG points with a projection value between 40 and 50
52	Percent of OG points with a projection value between 50 and 60
53	Percent of OG points with a projection value between 60 and 70
54	Percent of OG points with a projection value between 70 and 80
55	Percent of OG points with a projection value between 80 and 90
56	Percent of OG points with a projection value between 90 and 100
57	Percent of OG points with a projection value between 100 and 110
58	Percent of OG points with a projection value between 110 and 120
59	Percent of OG points with a projection value of 120
60	Percent of N points with a projection value below 40
61	Percent of N points with a projection value between 40 and 50
62	Percent of N points with a projection value between 50 and 60
63	Percent of N points with a projection value between 60 and 70
64	Percent of N points with a projection value between 70 and 80
65	Percent of N points with a projection value between 80 and 90
66	Percent of N points with a projection value between 90 and 100
67	Percent of N points with a projection value between 100 and 110
68	Percent of N points with a projection value between 110 and 120
69	Percent of N points with a projection value of 120
70	Percent of NZ points with a projection value below 40
71	Percent of NZ points with a projection value between 40 and 50
72	Percent of NZ points with a projection value between 50 and 60
73	Percent of NZ points with a projection value between 60 and 70
74	Percent of NZ points with a projection value between 70 and 80
75	Percent of NZ points with a projection value between 80 and 90
76	Percent of NZ points with a projection value between 90 and 100

77	Percent of NZ points with a projection value between 100 and 110
78	Percent of NZ points with a projection value between 110 and 120
79	Percent of NZ points with a projection value of 120
80	Percent of DU points with a projection value below 40
81	Percent of DU points with a projection value between 40 and 50
82	Percent of DU points with a projection value between 50 and 60
83	Percent of DU points with a projection value between 60 and 70
84	Percent of DU points with a projection value between 70 and 80
85	Percent of DU points with a projection value between 80 and 90
86	Percent of DU points with a projection value between 90 and 100
87	Percent of DU points with a projection value between 100 and 110
88	Percent of DU points with a projection value between 110 and 120
89	Percent of DU points with a projection value of 120

The projection value is the number of regularly-spaced 8 Å-long vectors emitted from each cavity points intercepting the protein surface (maximum of 120)

#### Correspondence between projection value and buriedness

Projection value	Buriedness, %
< 40	<33.3
40-50	33.3 - 41.6
50-60	41.6 – 50.0
60-70	50.0 - 58.3
70-80	58.3 - 66.6
80-90	66.6 – 75.0
90-100	75 .0 - 83.3
100-110	83.3 - 91.6
110-120	91.6 - 99.9
120	100

#### - Generating cavity-based pharmacophores

The **IChem --pharm** generates cavity-based pharmacophores ready to be used in BIOVIA, LigandScout or Shaper2.

Input directory: **\$ICHEM\_DIR/test/Volsite/pharm**

Input files: **protein.mol2, ligand.mol2**

Output files: **Pharmacophore.chm, Pharmacophore.mol2, Pharmacophore.plp, Pharmacophore.pml**

---

**> IChem --pharm -pCA 0 -pCZ 1.56 -pO 0 -pN 0 -pOD1 0 -pNZ 0 -pZn 3.49 volsite protein.mol2 ligand.mol2**

---

Usage of -p options (-p CA, -pCZ, -pO, -pN, -pOD1, -pNZ, -pZn) enables to filter out pharmacophoric features by PLP interaction energies (Gehlhaar et al. Chem. Biol. 1995, 2, 317-324) according to a

given threshold. We recommend using the above-described thresholds to optimize the accuracy of cavity-based pharmacophores.

4 files are generated: three in various pharmacophoric formats (BIOVIA .chm, LigandScout .pml, Shaper2 .mol2). The .plp file summarizes the PLP interaction energy (in kcal/mol) of each feature.



**Some notes:**

- The ligand may be a simple atom with user-defined coordinates!
- In ligand-restricted mode, druggability is only estimated for non-truncated full cavities (3<sup>rd</sup> column = ALL)
- The druggability model is only valid for standard VolSite parameters (step=1.5, boxes=20, b = 55, n =5)! If you change these parameters, keep in mind that the predicted druggability value has no meaning
- For large and accessible cavities (e.g. protein-protein interfaces), standard settings will output small-sized cavities. The most reliable approach to treat such flat cavities is to reduce the value of the *b* parameter (e.g. b= 45) until you get a set of cavity points that fills the expected interface area.
- If you want to consider accessory molecules (co-factor, nucleic acids, water) as part of your protein, use the corresponding options (--dna, --cofactor, --solvent, --hydrogen; see all arguments page 5)
- The additional --OutExclu parameter adds exclusion spheres on cavity-based pharmacophores created with the --pharm option.

## 6. PDB processing (pdbconv)

**Ichem pdbconv** is a tool to parse and process PDB files, automatically detect bound ligands (HET code) and their cavity, and estimates their druggability. It is the protocol that we currently use to set-up the sc-PDB database of druggable protein-ligand complexes (<http://bioinfo-pharma.u-strasbg.fr/scPDB>). For more information on the parsing process, please see: <http://bioinfo-pharma.u-strasbg.fr/scPDB/ABOUT>

***pdbconv*** *protein[.pdb|.mol2]* *output\_dir* *pdb\_id*

--wMOL2            Use MOL2 File as Input. PDB Options are not available  
 --wUnDrug        Output undruggable cavities  
 --noLig          PDB with no Ligand

By default all the following options are included. All chains will be kept

[PDB Options]

--HARMSIZE       Harmonize size line to 80 characters  
 --MSEMET        Change MSE to MET  
 --CSECYS        Change CSE to CYS  
 --MOVHET        move HETATM to the end of file  
 --ALTATM        select alternative atoms  
 --NUMATM        renumber atoms  
 --UPDMAS        update the MASTER line  
 --TOMOL2        convert to a molecular representation (instead of flat file)

if you use one of the option below, you MUST use also --TOMOL2 option or use --wMOL2 option

[MOL2 Options]

--RESTYP        apply Residue Class (cofactor/STD\_AA/MOD\_AA/Ligand ...)  
 --BONDSE        create bonds  
 --CLNUNW        clean unwanted residues  
 --MOL2TY        apply MOL2 types according to templates  
 --SPLITM        split molecule into protein/ligand/solvent  
 --SelChain N    List of chains to keep, separated by underscore  
 --SELWAT        select water molecules  
 --SELLIG        select ligand

**Input directory:** *\$ICHEM\_DIR/test/PDBCONV*

**Input file:** *2RH1.pdb*

**Output directory:** *./output*

---

**> IChem *pdbconv* *2RH1.pdb* *output* *2rh1***

---

For each detected ligand, a mol2 file is given along with the corresponding ligand-free protein (mol2 file), and Volsite cavities (mol2 file) in the **output** directory (to be created before executing the command). In addition, cavity descriptors are stored in the **\_descriptor.txt**.



## 7. Buried Surface Area calculation

The **IChem** *utils bsa* command computes the buried surface area of the ligand

*utils bsa protein ligand*

Input directory: `$ICHEM_DIR/test/BSA`

Input files: `protein.mol2`, `ligand.mol2`

Output file: `ligand.bsa`

---

**> IChem *utils bsa protein.mol2 ligand.mol2* > ligand.bsa**

---

The output is a 3 columns table:

REF	70.1934	347.125
-----	---------	---------

1st column: protein-ligand header

2<sup>nd</sup> column: buried surface area (%)

3<sup>rd</sup> column: ligand volume (Å<sup>3</sup>)

## 8. Protein-bound ligand fragmentation

The **IChem *utils frag*** command fragments in 3-D space a protein-bound ligand structure using a method aimed at detecting substituted ring cores. First, a ring perception algorithm is used to automatically detect aromatic and aliphatic rings. Acyclic atoms are then parsed to assign either a linker or substituent label, as whether to the corresponding bonds are connecting two rings or not. Linker atoms are left unchanged. In case of substituent atoms, single bonds involving the closest apolar carbon (in terms of bond distance) to any ring are later cleaved at the condition that the cleaved bond is at least 3 bonds away from the cyclic root atom. An anchoring atom (Z label) is then added to each of the remaining fragments to indicate the cleavage site

For more information, see:



Desaphy J, Rognan D.  
J Chem Inf Model. 2014 Jul 28;54(7):1908-18  
[sc-PDB-Frag: a database of protein-ligand interaction patterns for Bioisosteric replacements.](#)

***utils frag protein ligand***

**Input directory:** \$ICHEM\_DIR/test/Frag

**Input file:** protein.mol2, ligand.mol2

**Output files:** 2rh1\_CAU\_1\_protein\_FRAG\_1\_MOLE.mol2, 2rh1\_CAU\_1\_protein\_FRAG\_1\_INTS.mol2, xx.mol2

---

**> IChem *utils frag protein.mol2 ligand.mol2***

---

4 fragments are detected and their interaction with the target recorded in the terminal output. To be saved, a fragment should exhibit at least 5 interactions with the target.diff

Three mol2 files are outputted:

**2rh1\_CAU\_1\_protein\_FRAG\_1\_MOLE.mol2:** fragment from the original ligand

**2rh1\_CAU\_1\_protein\_FRAG\_1\_INTS.mol2:** interaction pseudoatoms (merged mode) for the fragment

**xxx.mol2:** interaction pseudoatoms (merged mode) for the entire ligand

### Some notes:



- Protein and ligand 3D coordinate should be in the same coordinates frame in mol2 file format.
- Several fragments can be generated. Each fragment has an index (header\_FRAG\_index.mol2) along with the corresponding interaction pseudoatoms (header\_FRAG\_index\_INTS.mol2)

## 9. Fingerprint similarity measures

The **IChem** **sims** command enables you to compute similarities with various metrics between two fingerprints.

### Fingerprint Similarity

**sims ref comp** (1)  
**sims file** (2)  
**sims RefInt Complnt** (3)

#### [General options]

--wlnts To use interactions instead of fingerprints (option 3 only)  
 --small Use small fingerprint (option 3 only)  
 --binary To add when the fingerprint is binary  
 -metric N (TC) Select the metric  
     TC Tanimoto metric (default)  
     HM Hamming distance  
     RT Ref Tversky  
     FT Fit Tversky  
     DI Dice  
     SO Soergel

Input directory: **\$ICHEM\_DIR/test/SIMILARITY**

Input file: **FP1.txt, FP2.txt, FP.txt**

Output files: **sim.txt, matrix.txt**

---

**> IChem --binary sims FP1.txt FP2.txt >sim.txt**

---

Fingerprints of diverse nature are accepted

STD	Standard (1 0 21 0 0 3)
SVM	SVM format (1:1 3:21 6:3)
CMP	Compressed (1 [1 21 [2 3)

The output is a table listing the reference fingerprint, the target fingerprint and the similarity score:

FP1	FP2	0.809524
-----	-----	----------

Alternatively, a full similarity matrix (all against all) can be generated by concatenating all fingerprints to compare (see FP.txt input file) and running IChem in matrix mode

---

**> IChem --binary sims FP.txt >matrix.txt**

---

The output is a table listing the reference fingerprint, the target fingerprint and the similarity score as follows:

FP1	FP1	1
FP1	FP2	0.809524
FP1	FP3	0.652174
FP2	FP1	0.809524
FP2	FP2	1
FP2	FP3	0.714286
FP3	FP1	0.652174
FP3	FP2	0.714286
FP3	FP3	1

## 9. Detection and analysis of protein-protein interfaces (DetectPPI)

DetectPPI is a tool to detect and analyze protein-protein interfaces, and outputs cavities remote to every detected PPI

For more information, see:



Da Silva, F., Desaphy, J., Bret, G. and Rognan, D. (2015) [IChemPIC: A Random Forest Classifier of Biological and Crystallographic Protein-Protein Interfaces](#). J. Chem. Inf. Model, 55, 2005–2014.

**Input directory:** \$ICHEM\_DIR/test/PPI

**Input file:** 4NN6.pdb

**Output files:**

4NN6_interaction_A_B.ints	4NN6_prot_B_C.mol2	CAVITY_all_N2_ALL.mol2	CAVITY_B_N1_ALL.mol2
4NN6_interaction_A_C.ints	4NN6_site_A_B.mol2	CAVITY_all_N3_ALL.mol2	CAVITY_B_N2_ALL.mol2
4NN6_interaction_B_C.ints	4NN6_site_A_C.mol2	CAVITY_all_N4_ALL.mol2	CAVITY_B_N3_ALL.mol2
4NN6_ints_A_B.mol2	4NN6_site_B_C.mol2	CAVITY_all_N5_ALL.mol2	CAVITY_C_N1_ALL.mol2
4NN6_ints_A_C.mol2	CAVITY_all_N10_ALL.mol2	CAVITY_all_N6_ALL.mol2	descriptor.sre
4NN6_ints_B_C.mol2	CAVITY_all_N11_ALL.mol2	CAVITY_all_N7_ALL.mol2	VolSite_Stat.csv

---

> **IChem detectppi 4NN6 4NN6.pdb**

---

IChem requires a PDB file with orientated polar hydrogen atoms. You can use any program for that specific task although we strongly recommend the usage of ProToss (<http://protoss.zbh.uni-hamburg.de/>) for that purpose.

PDB\_interaction\_X\_Y.ints: interactions between protein chains X and Y

PDB\_ints\_X\_Y.mol2: interaction pseudoatoms (centered mode) between protein chains X and Y

PDB\_prot\_X\_Y.mol2: mol2file of protein chains X and Y

PDB\_site\_X\_Y.mol2: mol2file of interacting residues from chains X and Y

CAVITY\_all\_X\_ALL.mol2: Volsite cavities for the entire protein assembly

CAVITY\_X\_NY\_ALL.mol2: Volsite cavities in chain X only

descriptor.sre: set of 45 descriptors of each PPI

VolSite\_Stat.csv: Volsite output

Here is the meaning of the 45 PPI descriptors:

Name	Description
nPTS	Total number of interaction points
Hydro	% of hydrophobic interaction points
Aro	% of aromatic interaction points
Hbond	% of hydrogen-bond interaction points
Ionic	% of ionic bond interaction points
Hydro1	% of hydrophobic points (25 %<Burial <33.3%)
Hydro2	% of hydrophobic points (33.3 %<Burial <41.6%)
Hydro3	% of hydrophobic points (41.6 %<Burial <50%)
Hydro4	% of hydrophobic points (50 %<Burial <58.3%)
Hydro5	% of hydrophobic points (58.3 %<Burial <66.6%)
Hydro6	% of hydrophobic points (66.6 %<Burial <75%)
Hydro7	% of hydrophobic points 75 %<Burial <83.3%)

<b>Hydro8</b>	% of hydrophobic points (83.3 %<Burial <91.6%)
<b>Hydro9</b>	% of hydrophobic points (91.6%<Burial <100%)
<b>Hydro10</b>	% of hydrophobic points (Burial =100%)
<b>Aro1</b>	% of aromatic points (25 %<Burial <33.3%)
<b>Aro2</b>	% of aromatic points (33.3 %<Burial <41.6%)
<b>Aro3</b>	% of aromatic points (41.6 %<Burial <50%)
<b>Aro4</b>	% of aromatic points (50 %<Burial <58.3%)
<b>Aro5</b>	% of aromatic points (58.3 %<Burial <66.6%)
<b>Aro6</b>	% of aromatic points (66.6 %<Burial <75%)
<b>Aro7</b>	% of aromatic points 75 %<Burial <83.3%)
<b>Aro8</b>	% of aromatic points (83.3 %<Burial <91.6%)
<b>Aro9</b>	% of aromatic points (91.6%<Burial <100%)
<b>Aro10</b>	% of aromatic points (Burial =100%)
<b>Hbond1</b>	% of hydrogen bond points (25 %<Burial <33.3%)
<b>Hbond2</b>	% of hydrogen bond points (33.3 %<Burial <41.6%)
<b>Hbond3</b>	% of hydrogen bond points (41.6 %<Burial <50%)
<b>Hbond4</b>	% of hydrogen bond points (50 %<Burial <58.3%)
<b>Hbond5</b>	% of hydrogen bond points (58.3 %<Burial <66.6%)
<b>Hbond6</b>	% of hydrogen bond points (66.6 %<Burial <75%)
<b>Hbond7</b>	% of hydrogen bond points 75 %<Burial <83.3%)
<b>Hbond8</b>	% of hydrogen bond points (83.3 %<Burial <91.6%)
<b>Hbond9</b>	% of hydrogen bond points (91.6%<Burial <100%)
<b>Hbond10</b>	% of hydrogen bond points (Burial =100%)
<b>Ionic1</b>	% of ionic bond points (25 %<Burial <33.3%)
<b>Ionic2</b>	% of ionic bond points (33.3 %<Burial <41.6%)
<b>Ionic3</b>	% of ionic bond points (41.6 %<Burial <50%)
<b>Ionic4</b>	% of ionic bond points (50 %<Burial <58.3%)
<b>Ionic5</b>	% of ionic bond points (58.3 %<Burial <66.6%)
<b>pwlIonic6</b>	% of ionic bond points (66.6 %<Burial <75%)
<b>Ionic7</b>	% of ionic bond points 75 %<Burial <83.3%)
<b>Ionic8</b>	% of ionic bond points (83.3 %<Burial <91.6%)
<b>Ionic9</b>	% of ionic bond points (91.6%<Burial <100%)
<b>Ionic10</b>	% of ionic bond points (Burial = 100%)

Please note that we recommend the usage of our web interface (<http://bioinfo-pharma.u-strasbg.fr/IChemPIC>) to predict the relevance (biological, crystallographic) of any possible interface from a PDB file.